

SWISS LITHOPUMP®

→ FOR ALL LITHOCLAST® DEVICES



EMS+

**VENTURI
PRINCIPLE**

PRECISELY ADJUSTABLE → FOR ALL SUCTION REQUIREMENTS OF LITHOTRIPSY PROCEDURES

FAILSAFE AND ULTRAQUIET

- > Runs on compressed air (venturi principle, no electricity)

EASE OF USE

- > Plugs into the compressed air connection of the Swiss LithoClas[®]

CLEAN AND HYGENIC

- > Double bacteria retention filter for increased hygiene safety
- > Minimum air consumption due to pneumatic auto control circuit

SWISS LITHOPUMP[®]

Swiss LithoPump[®] basic unit including compressed air connection, filters and connection tubing

FT-184#A

- > Swiss LithoPump[®] with Swiss LithoClas[®] Master



ACCESSORIES

Overflow protection filters (10x) **DT-096**

Disposable suction bags 2l (25x) **DT-097**

Double reusable suction jar set 2l **EP-035**

Autoclavable suction jar **EP-036**

Reusable suction jar **EP-037**

(for disposable suction bags)

Clamp holder **EP-038**

Change-over valve **EG-088**

Cart **DW-031**



EMS Electro Medical Systems S.A.
Chemin de la Vuarpillièrre 31
1260 Nyon
Switzerland

Tel. +41 22 99 44 700
Fax +41 22 99 44 701
welcome@ems-ch.com
www.ems-urology.com

EMS 