

SWISS LITHOPUMP®

→ DEDICATED SUCTION SOLUTION
FOR SWISS LITHOCLAST®



EMS 

**VENTURI
PRINCIPLE**

PRECISELY ADJUSTABLE → FOR ALL SUCTION REQUIREMENTS OF LITHOTRIPSY PROCEDURES

FAILSAFE AND ULTRAQUIET

- > Runs on compressed air (venturi principle, no electricity)

EASE OF USE

- > Plugs into the compressed air connection of the Swiss LithoClast®

CLEAN AND HYGENIC

- > Double bacteria retention filter for increased hygiene safety
- > Minimum air consumption due to pneumatic auto control circuit

SWISS LITHOPUMP®

Swiss LithoPump® basic unit including compressed air connection, filters and connection tubing

FT-184#A

- > Swiss LithoPump® with Swiss LithoClast® Master



ACCESSORIES

Overflow protection filters (10x)

DT-096

Disposable suction bags 2l (25x)

DT-097

Double reusable suction jar set 2l

EP-035

Autoclavable suction jar

EP-036

Reusable suction jar

EP-037

(for disposable suction bags)

Clamp holder

EP-038

Change-over valve

EG-088

Cart

DW-031



EMS Electro Medical Systems S.A.
Chemin de la Vuarpillière 31
1260 Nyon
Switzerland

Tel. +41 22 99 44 700
Fax +41 22 99 44 701
welcome@ems-ch.com
www.ems-urology.com

EMS 